


## Overview

### Offline



# Changing in offline analysis

## SSVEP\_ML.ipynb

1. Changing Filter Method

```
from scipy import signal

def butter_bandpass(lowcut,highcut,fs,order):
    nyq = 0.5*fs
    low = lowcut/nyq
    high = highcut/nyq
    b,a = signal.butter(order,[low,high],'bandpass')
    return b,a

def butter_bandpass_filter(data,lowcut = 6, highcut = 30, order = 4, axis = 1):
    b,a = butter_bandpass(lowcut,highcut,512,order)
    y = signal.filtfilt(b,a,data,axis=axis)
    return y

Epochs_data = butter_bandpass_filter(Epochs.get_data(), lowcut = 2, highcut= 40, axis = 2)
```

### 2. Saving ML model

```
y_train = train_label # Get true label
y_test = test_label

svm_model = SVC(C = 1, kernel= 'linear') # Using a linear kernel
svm_model.fit(x_train, y_train)

print(x_train.shape)

print('accuracy', svm_model.score(x_train, y_train))
label_names = ['12Hz', '6Hz', '24Hz', '30Hz']

with open("trained_model/SVM_model.pkl", "wb") as file:
    pickle.dump(svm_model, file)

GetConfusionMatrix(svm_model, x_train, x_test, y_train, y_test, label_names)
```

## **EEG** data Streaming

### biosemi\_stream.ipynb

1. Select target EEG channels

```
raw_eeg = mne.io.read_raw_fif("datasets/biosemi_SSVEP.fif")
select_ch = ['01','0z','P03','P04','P0z','Pz']
raw_eeg = raw_eeg.pick(select_ch)
```

2. Change name of StreamInfo and chunk\_size same as your sampling rate

## **Real-time BCI-UI**


## Main.py

Config parameters

Don't forget to match your stream name

```
ch_names = ['01','0z','P03','P04','P0z','Pz']
num_channels = len(ch_names)
samp freq = 512
window size second = 4
if __name__ == '__main__':
   mp.freeze support() # Ensure compatibility with multiprocessing on Windows
   # Use multiprocessing queues for status updates
   status queue = mp.Queue()
   queue1 = mp.Queue()
   command queue = mp.Queue()
   # Initialize the Model and View with LSL streaming
   model = EEGModel(num channels=num channels, samp freq = samp freq, window size second = window size second, band pass = (2,40),
                    stream name="SomSom", status queue=status queue, command queue=command queue)
   view = RealTimeView(model, ch_names, samp_freq=samp_freq, window_size_second=window_size_second)
   # Start the streaming process in the model
   model.start_streaming()
   # # Run the DearPyGUI rendering loop in the main thread
   view.setup windows() # Setup both windows
   view.render_loop()
                         # Start the rendering loop
   # Stop the model when exiting
   model.stop streaming()
```

### View.py



# **Unity**

If you want to receive command in Unity, copy this code to your C# script

#### Print Log in callback function

```
using UnityEngine;
using UnityEngine.UI;
using System;
using System.Text;
using System.Net;
using System.Net.Sockets;
using System.Threading;
public class PythonToUnity1 : MonoBehaviour
   Thread receiveThread;
   UdpClient client;
   public int port = 1880;
   private bool startRecieving = true;
   private string data;
   // Start is called before the first frame update
   void Start()
       receiveThread = new Thread(
            new ThreadStart(ReceiveData));
       receiveThread.IsBackground = true;
       receiveThread.Start();
   public void ReceiveData()
       client = new UdpClient(port);
       while (startRecieving)
            try
                IPEndPoint anyIP = new IPEndPoint(IPAddress.Any, 0);
                byte[] dataByte = client.Receive(ref anyIP);
                data = Encoding.UTF8.GetString(dataByte);
                Debug.Log(data);
            catch (Exception err)
                print(err.ToString());
   // Update is called once per frame
   void Update()
```

# Files Running Steps

- 1. Run biosemi\_stream.ipynb
- 2. Run Main.py

